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INFORMATION MEDICALE IMPORTANTE  
EXRAY<sup>®</sup>20 EXRAY<sup>®</sup>30 et EXRAY<sup>®</sup>35

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Reference: SUP\_7.016 version 11  
Date of last update: 01/2024

Date of first EC marking:  
-EXRAY<sup>®</sup> 20: 2012  
-EXRAY<sup>®</sup>30 / EXRAY<sup>®</sup>35: 2012

**MDSS-UK RP Ltd.**  
8 Wilmshurst Road  
Rushmore  
Manchester M14 5TP  
UNITED KINGDOM

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**CAUTION: USA Federal law restricts this device to sale by or on the order of physician.**

## FDA Specific - CAUTION

**Specific IFU is dedicated to USA.**  
Please refer to our website, section "IFU": [www.synchromedical.com](http://www.synchromedical.com).  
In case you prefer receiving IFU in a paper-based version, please: send an email to [office@synchromedical.com](mailto:office@synchromedical.com) or call +33 428 710 320. IFU will be sent to the address you will provide without any additional fee.

## PURPOSE

EXRAY<sup>®</sup> screws are designed to relieve pain and disability of the forefoot by attaching and stabilising bone segments in elective osteotomies on metatarsal bones and phalanges of the foot, in order to optimise the correct bone fusion of the segments concerned. These devices are only intended for use for the forefoot of a mature skeleton.

## DESCRIPTION:

The EXRAY<sup>®</sup> screw system is intended for single use and is comprised of EXRAY<sup>®</sup> 20, EXRAY<sup>®</sup> 30 and EXRAY<sup>®</sup> 35 screws.

EXRAY<sup>®</sup> 20 screws comprise a single part implant in the form of a pin fitted with a threaded portion and a connector portion with a narrower diameter. EXRAY<sup>®</sup> 20 screws are non-compressive 2mm screws with a continuous thread.

EXRAY<sup>®</sup> 30 screws comprise a single part implant. These 3mm self-compressive headless screws are available in several different lengths.

EXRAY<sup>®</sup> 35 screws comprise a single part implant. These 3.5mm self-compressive headless screws are available in several different lengths.

EXRAY<sup>®</sup> screws are manufactured in PEEK according to Standard ASTM F2026.

## INDICATIONS, CONTRAINDICATIONS AND POSSIBLE SIDE EFFECTS

**Indications:**  
The EXRAY<sup>®</sup> screw system is indicated as an attachment system for small bone fractures or minor bone reconstruction in mature skeleton in the following cases:

- noncoroncal or bicoronal osteotomies of the foot
- metatarsal or distal or proximal metatarsal osteotomies
- fusion of the first metatarsophalangeal and interphalangeal joint - attachment in the case of treatment of Hallux Valgus (e.g. scarf osteotomy, chevron osteotomy, etc.)
- Akin osteotomy

**CAUTION: to be used by or on the order of a surgeon.** The surgeon must take note of the documents accompanying the device. No specific training is required for the understanding and use of the device. The surgeon's qualifications and the reading of the accompanying documents are enough.

## Contraindications:

- A non-exhaustive list of contraindications is as follows:
- Any sign of generalised or local infection.
  - Pathological obesity.
  - Pregnancy.
  - Any other medical or surgical condition that may compromise the success of surgery with instruments, such as the presence of malignant tumours, or serious congenital anomalies, an increase in sedimentation rates that cannot be attributed to other diseases, an increase in the number of white blood cells or a downward trend in such blood cells.
  - Suspected or known allergy or intolerance to the implant's component materials.
  - Any situation requiring the use of different materials
  - Any case not listed in the indications.
  - Any patient who is not willing to follow the postoperative instructions.
  - Any patient in whom use of the implant may interfere with anatomical organs or some expected physiological function. The contraindications related to these devices are similar to those related to other osteosynthesis instruments. These osteosynthesis instruments have not been designed for, intended or sold for any use other than those indicated.

## Possible side effects:

- In addition to risks associated with any surgery to the extremities, the list of possible undesirable side effects is as follows (non-exhaustive list):
- Migration, fracture / rupture, expulsion, dislocation, or undesirable movement of one or more of the components
  - Persistent pain and disability
  - Bone loss related to short-circuiting of the constraints.
  - The body's reaction to a foreign object owing to the presence of the implant (e.g. the development of a tumour, auto-immune disease and / or poor healing, inflammation and/or irritation)
  - Skin incompatibility with components where there is insufficient tissue coverage with pain and abnormal sensations due to the volume of the material.
  - Superficial or deep infection with inflammatory reaction.
  - Septis
  - Pseudarthrosis or delay of bone fusion.
  - Recurrence of the deformation, loss of correction
  - Excessive loss of blood during the operation and / or haematoma.
  - Lesions and / or vascular disorders (thrombosis)
  - Peripheral circulatory complication (thrombosis) and / or pulmonary embolism.
  - Inability to resume the activities of normal everyday life.
  - Compression of nearby tissues or organs.
- NOTE: Some additional surgical procedure may be required to remedy some of these possible undesirable events.

## INSERTION OF THE DEVICE

The EXRAY<sup>®</sup> screw system must be used with instruments specially designed for the device according to the surgical technique concerned. It must not be used with any other instruments, unless this has specifically been recommended in other ADSM documentation, as using the system with other instruments poses the risk of incompatibility and cannot be guaranteed. ADSM shall not be held responsible in the case of using instruments that are not related to the implanted system. NEVER REUSE AN IMPLANT. Reuse of an implant presents the risk of contamination and loss of mechanical integrity.

## PACKAGING

The packaging of the device must be intact on receipt. Any device that has been damaged or any device whose packaging is damaged must not be used and must rather be returned to ADSM.

## STORAGE AND HANDLING CONDITIONS

The implants must be stored in their original packaging, out of direct sunlight, humidity and extreme temperatures.

## MAGNETIC RESONANCE

Magnetic Resonance Imaging (MRI) is a medical imaging technique that provides non-invasive views of the body interior using a strong magnetic field. Consequently, any ferromagnetic metallic foreign body is to be proscribed which can cause burns or injuries by displacement of the body under the effect of attraction of the magnetic field. The PEEK used to manufacture the EXRAY<sup>®</sup> screw system is therefore a non-ferromagnetic and non-conductive polymer (plastic). Consequently, it poses no known risk of possible heating or migration due to exposure to the magnetic resonance environment.

## VERIFICATION

The devices must always be checked before use. Those that present signs of damage or scratches on the surface must not be used. EXRAY<sup>®</sup> implants are delivered sterile and have been sterilised through gamma radiation. You need to carefully check that the round disc (sterilisation mark) on the packaging is red. Do not use an implant whose disc isn't red. You also need to check that the expiry date on the packaging has not passed. ADSM disengages itself from any liability in the case of the use of implants past their expiry date. Any implant with damaged or torn packaging must not be resterilised but rather be returned to ADSM. TAKE NOTE: Do not resterilise an implant that has already been sterilised.

## WARNINGS AND PRECAUTIONS

Every surgery carries the risk of failure, as many external factors can compromise the outcome, especially if the patient's state of health is already compromised. The implants must only be used in cases where bone fusion is deemed to be possible. No implant has been designed to withstand body weight in the long term without being accompanied by solid bone fusion. If a delay or lack of consolidation occurs, the implant may migrate and/or rupture. The functional product lifetime is of 6 months. In fact, the first bone fusion occurs around 90 days and keeps progressing. Therefore, from 6 months the product is still implanted but it has no functional performance. Given the current state of knowledge and device anteriority, the risk caused by aging cannot be totally ruled out. The following of preoperative and operative procedures, knowledge of surgical techniques, the selection of implants and their correct placement as well as appropriate selection and the patient's cooperation are factors leading to the operation's success. It has been shown that patients who smoke present a high risk of failure to obtain bone fusion. They need to be informed of this and any possible consequences. Obese patients, those who are malnourished and/or consume alcohol as well as patients with poor muscle or bone quality and/or those suffering from paralysis are also poor candidates for osteosynthesis surgery.

## OTHER PREOPERATIVE, PEROPERATIVE AND POSTOPERATIVE WARNINGS

**Implant selection**  
Selecting the suitable size of the implant for each patient is a crucial factor in the success of the operation, and it is important to refer to the surgical technique in this regard. Once implanted, the implant is subjected to repeated strain, and its strength is limited

by adaptation of its geometrical structure to the size and shape of human bone. In order to minimise this strain, it is important to pay careful attention to the patient's selection criteria, the correct placement of the implant, and postoperative care.

In fact, conversely, these constraints can lead to excessive strain placed on the material, resulting in deformation, rupture or loosening of the device likely to cause damage or the need to remove the implant prematurely. Any use in body areas not recommended by ADSM cannot be guaranteed.

## Preoperative precautions

- Only patients who meet the criteria described in the indications should be selected.
- Patients who meet the criteria described in the contraindications should not be selected.
- The implants must be handled and stored with the utmost care. They should not be scratched or damaged.
- Before carrying out the procedure, the surgeon must check availability of the range of implants and associated instruments. They must personally handle all the devices before use in order to familiarise themselves with insertion of the implant.
- Additional sterile components must be available in the event of unexpected need.

## Preoperative precautions

- Instructions concerning the surgical technique should be followed closely. Any breakage or improper use of the implant may injure the patient or surgical team.
- Specific instrumentation must be used to insert the device.
- Take note: Do not use a device of the incorrect size (diameter and/or length), as this may damage the surrounding tissue.
- Where possible and necessary, an imaging system may be used to check on placement of the device.

## Preoperative precautions

- The instructions and postoperative warnings provided to the patient by the doctor, and the patient's ability to follow them are extremely important, so as not to compromise the desired bone fusion.
- Detailed instructions on the device's limits need to be provided to the patient. The patient must be informed that any deformity, loosening, and/or rupture of the device are complications that may occur during the postoperative rehabilitation phase following early use or excessive muscle activity. The risk of the onset of these complications could be heightened if the patient is active, frail, mentally ill or incapable of using some immobilisation device or external support. The patient must be advised in order to prevent falls, jumps or hard impacts on the area operated on or about any other movement that could compromise the desired bone fusion.
  - To increase the likelihood of success of the surgery, the patient must not expose themselves to excessive loads which could cause weakening or rupture of the device. The patient must be informed of this risk, and advised so they are aware of how to limit their physical activity, in particular twisting movements as well as any other participation in some sporting activity.
  - The patient should be informed about avoiding smoking or consuming alcohol during the bone consolidation process.
  - Lack of bone fusion would lead to excessive, repeated constraints placed on the implant. Through the mechanism of fatigue, these constraints could end up causing some deformation, loosening or rupturing of the device. It is important to immobilise the fusion area and to check consolidation through x-ray. Should lack of bone fusion persist, or in the case of rupture, the device will need to be adjusted and/or removed immediately, before any serious injury occurs.
  - Any implant collected following ablation must be treated in a manner to prevent its reuse in some other surgical procedure in accordance with the regulations described in the removal and disposal section for medical devices contained in this manual.

## ADDITIONAL INFORMATION

For any additional information about the device or to enquire about specific surgical techniques, kindly contact ADSM's client services or the distributor.

The instruction for use can be downloaded from the [www.adsm-group.com/ifu](http://www.adsm-group.com/ifu). If the website is unavailable or you face difficulties to access the document, ADSM can send the instruction for use either by e-mail in PDF format or by post at no extra cost (see contact details above for sending us your request).

## PRODUCT COMPLAINTS

Any customer or user of this system who has a complaint or any reason for dissatisfaction concerning the quality of a product, its identity, durability, reliability, safety, effectiveness and/or its performance must notify the distributor or ADSM. In addition, should any of the implanted components not function well, or if poor functioning is suspected, the distributor or ADSM must be notified about this. If one of ADSM's products is ever suspected of having caused or contributed towards a patient's death or serious injury, ADSM or the distributor and the competent authority of the Member State in which the user and/or patient is established should be informed of this immediately by telephone, fax or written correspondence. For any complaint, please indicate the name, reference number and batch number of the component(s) concerned, your name and address, the nature of your complaint, and specify whether a written report is requested from the distributor or ADSM.

## INFORMATION TO BE PROVIDED TO THE PATIENT

The EXRAY<sup>®</sup> screw system has been designed to attach and stabilise bone segments in elective osteotomies on metatarsal bones and phalanges of the foot, in order to optimise the correct bone fusion of the segments concerned. The device has only been designed for use in the forefoot. Specific instructions for the implant are described in the paragraph: "INDICATIONS, CONTRAINDICATIONS AND POSSIBLE UNDESIRABLE EFFECTS". This device does not restore the functionality that can be expected of healthy bone and the patient should not expect unrealistic functional results. In addition, the anatomy of the human body limits the size of any artificial consolidation device used in surgery. This geometrical limitation increases the possibility of the onset of mechanical complications such as dismantling, deformation or rupturing of

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