-FR-

RECOMMANDATIONS

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SynchroMedical
The international reference

7b rue Lavoisier
69680 Chassieu - FRANCE
Tél : +33 (0)4 28 71 03 10
Fax : +33 (0)4 28 71 03 20
office@synchronmedical.com

UK CA 

UK-RP : MDSS-UK RP Ltd.
6 Wilmslow Road, Rusholme
Manchester M14 5TP
United Kingdom

CH REP MDSS CH GmbH
Laurenzenvorstadt 61
5000 Aarau
Switzerland

-EN-
RECOMMENDATIONS

***FOR CLEANING AND STERILISING
INSTRUMENTS AND ANCILLARIES***

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Introduction

This document has been developed to provide recommendations on the cleaning and sterilization of reusable ADSM instruments for multiple patients.

The instructions supplied in this notice apply to all ADSM instruments and ancillaries.

It is incumbent on a health care centre to ensure that the re-treatment is done using appropriate equipment and materials and that the personnel responsible for the re-treatment have been correctly trained in order to obtain the desired result. The equipment and the processes must be validated and regularly monitored.

Warnings

- Surgical instrumentation must be manipulated and/or used by trained, qualified persons having read the present instruction notice
- Surgical instrumentation must be used in premises compatible with standard aseptic conditions and with current practices in health care establishments
- The ancillaries and instruments are delivered clean but NON-STERILE.
- Single use devices must not be re-sterilised with the exception of implants delivered non-sterile. These instructions do not apply to ADSM medical devices supplied sterile which must not be reused
- ADSM ancillaries and instruments are designed exclusively for the implantation or the explantation of ADSM implants. Under no circumstances can they be used for surgical instruments from other manufacturers
- The transfer of an ancillary between two health care establishments is forbidden without prior return to ADSM
- Certain ancillaries are constituted by temporary invasive devices such as drillbits, reamers, millbits, taps and probes. These instruments must be used with precaution.
- The equipment must only be used for the function for which they are dedicated in operating techniques. Any deterioration of the equipment can lead to risks of malfunction

Generalities

- Before being used, the instruments, whether new or not, must be treated carefully according to the instructions
- The instruments boxes do not constitute a sterile barrier; a sterile packaging must be used to keep them in a sterile state
- Ancillary equipment is delivered non-sterile in transport containers. A transfer sheet (which specifies the non-sterile state of the ancillary) is supplied with each ancillary for checking at reception and then before sterilisation
- Instruments or ancillaries generally have a long service life: however, poor handling or unsuitable protection can rapidly reduce it. Instruments which are no longer efficient because of long use, poor handling or inappropriate maintenance, must be returned to ADSM
- Any malfunction concerning the instruments must be indicated to ADSM
- ADSM instruments and ancillaries must undergo all the steps of decontamination, cleaning, disinfection, inspection and final sterilisation before being sent back to ADSM. The documentation on decontamination must be supplied with instruments returned to ADSM
- Operating techniques for each of the families of implants of the ADSM range are available in order to provide the surgeon with further information concerning the surgical technique. This documentation must be read
- Any non-sterile medical device received must be inspected upon receipt to verify the integrity and the proper functioning of the medical device.

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- Washer-disinfectors used by the health care facility must comply with ISO 15883.

Cleaning

Caution

- The following precautions must be observed by hospital personnel who work with contaminated or potentially contaminated medical devices. Pointed or cutting devices must be handled with great care
- Do not use metal brushes or scouring pads during manual cleaning procedures. These items risk damaging the surface and the finish of the instruments. Use nylon brushes with soft bristles and swabs with dimensions adapted to the devices to treat
- Do not allow contaminated devices to dry before re-treatment so as to facilitate later steps
- Saline solution and cleaning/disinfecting agents containing aldehyde, mercury, active chlorine, chloride, bromine, bromide, iodine or iodide are corrosive and must not be used
- Do not use mineral oil or silicone lubricants as they coat the micro-organisms preventing direct contact of the steam with the surface and are difficult to eliminate.
- Cleaning and neutral pH enzymatic agents are recommended and preferred for cleaning reusable ADSM devices
- Only agents with a neutral pH must be used for cleaning aluminium devices. Contact with strong alkaline detergents or with solutions containing caustic soda, iodine or chlorine must be avoided as aluminium and stainless steel can be chemically attacked and the device damaged. The user must always refer to and comply with the instructions provided by the cleaning product manufacturer
- The use of hard water must be avoided. Softened tap water can be used for initial rinsing. Final rinsing must be done with purified water in order to eliminate deposits on the instruments. One or several of the following processes can be employed to purify water: ultrafiltration (UF), reverse osmosis (RO); de-ionisation or equivalent

Comment 1: Drills, broaches, rasps and cutting instruments must be carefully inspected after treatment with alkaline detergents in order to ensure that the cutting edges are not deteriorated.

Comment 2: In the case of first use of the instruments, the establishment shall perform a series of cleaning operations in order to eliminate any earlier contamination of the instruments. A minimum of 3 consecutive cleaning operations is recommended.

Instructions

Where the device will be used

- Remove excess liquid and organic tissues on the instruments with a lint-free disposable pad. Place the devices in a tray of distilled water and cover with a damp cloth
- Remove the plastic films and potential protective foam
- Visually check that all the instruments cited on the attached check-list are present

Comment: the instruments must be cleaned within 30 minutes following their use to minimise the possibility of drying before cleaning.

Confinement and transport

The instruments used must be sent to the procurement service in closed or covered containers in order to avoid any new risk of contamination

- When this is applicable, the instruments with several components must be disassembled for efficient cleaning. Make sure not to lose small screws and components
- When this is applicable, articulated instruments must be open

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- All cleaning agents must be prepared observing the dilution and the temperature recommended by the manufacturer. Softened tap water can be used for preparing the cleaning agents
- The use of recommended temperatures is important for optimal performances of the cleaning agents

Comment: new cleaning solutions must be prepared.

Procedure for automatic cleaning/disinfection

Soaking (pre-cleaning)

- Completely immerse the instruments in an alkaline solution having fungicide, bactericides and virucide powers and leave them to soak for at least 15 minutes at room temperature. We recommend the use of the product NEODISHER®SEPTO PRECLEAN ZP ([C] = 5 ml / L) of the brand DR WEIGERT Use a nylon, flexible bristle brush to softly brush the device until all visible stains are eliminated. Pay particular attention to asperities, cavities, contact surfaces, connectors and other difficult to clean areas. Action the articulated parts in order to get to difficult access areas

Comment: it is not recommended to prolong the soaking step for an indeterminate period at the risk of deteriorating the instruments.

- Remove the device from the alkaline solution and rinse abundantly with lukewarm water of the network for at least 5 minutes. Carefully and effectively rinse cavities, holes and any other difficult access areas

Cleaning

- Place the instruments in the tray of a suitable washing machine/autoclave and treat using the standard cycle of the instrument washing machine/autoclave as per the instructions from the manufacturer of the washing machine/autoclave
- Washing with a detergent and disinfectant solution suitable for the washing machine/autoclave (observe the concentrations and the temperatures of use of the products). We recommend the use of the product NEODISHER®SEPTO CLEAN ([C]=5 ml/L) DR WEIGERT's brand.

Disinfection

- Thermal disinfection assuring a value of A_0 greater than 3000, i.e. A minimum of 3 min at 93°C (or any other equivalent time/temperature pair)
- Rinse with demineralised water
- Drying not exceeding 110°C. We recommend the use of MEDIKLAR® SPECIAL drying activator ([C] = 0.8 ml / L) of the DR WEIGERT brand.

Comment: The ancillary equipment is designed to operate at room temperature. It is recommended not to use the equipment directly after steam sterilisation and to wait for the temperature to return to ambient. Before the intervention, it is the responsibility of the personnel of the operating theatre to check the proper coupling of the ADSM instruments with the equipment available in the operating theatre.

Inspection before packaging

- Carefully inspect each device to ensure the elimination of any visible contamination. If there is any contamination, repeat the cleaning/disinfection process
- Check the action of movable parts (e.g. hinges, connectors, sliding parts, springs, etc.) over their complete amplitude

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- Check that there is no deformation of the instruments, in particular instruments which are connected to motors
- When the instruments are assembled, check that the devices fit together correctly with any components in contact
- Check that there is no excessive play between the assembled components

Packaging

- Medical quality packaging or steam sterilisation sachets available commercially, can be used for wrapping individual instruments. The wrapping must be prepared using a protective double wrapping technique or an equivalent method
- For sterilisation, trays and boxes with covers can also be placed in an approved sterilisation container with sealing lid. Follow the sterilisation container manufacturer's instructions for insertion and the replacement of sterilisation filters in sterilisation containers
- Instrument transport and storage boxes and trays must be used in the following conditions:
 - All devices must be properly arranged in order to ensure contact of the steam with all instrument surfaces. The instruments must not be stacked or placed in close contact (use the arranging supports or silicone mat provided for this)
 - The user must check that the contents of the instrument box do not spill out, once the devices are arranged in the box
 - Only devices manufactured and/or distributed by ADSM must be placed in ADSM trays

Sterilization

- The cleaning/decontamination step is compulsory before sterilising. It is important that sufficient cleaning is done before sterilisation without which, the sterilisation can be compromised
- After cleaning, the instruments must be rapidly and meticulously reassembled and placed in their respective trays/locations in order to avoid any possible damage or contamination
- The instrument trays and supports must have been decontaminated and cleaned according to the same arrangements as those for the instruments
- The health care establishment is responsible for the internal procedures for reassembly, inspection and packaging of instruments in order to ensure the sterilising penetration of the steam and appropriate drying. The arrangements to make for the protection of sharp or potentially injuring edges of instruments must also be recommended by the health care establishment
- Sterilisation by steam/moist heat is the recommended method for the ADSM set of orthopaedic instruments. Refer to the table below for the parameters of the sterilisation cycles validated by ADSM ensuring a sterility assurance level (SAL) of 10^{-6}
- Sterilisation methods using ethylene oxide or plasma must not be used
- The recommendations of the manufacturer of the steriliser must always be followed. During the sterilisation of several instrument sets in a sterilisation cycle, ensure that the maximum load indicated by the manufacturer is not exceeded
- The set of instruments must be correctly prepared and packaged in the trays and/or boxes to allow the steam to penetrate and enter in direct contact with all surfaces

Validated parameters for steam sterilization:

	Type of cycle	Temperature (°C) / (°F)	Exposure Time	Drying Time
US	Pre-vacuum	132°C / 270°F	4 min	20 min
International*	Pre-vacuum	134°C / 273°F	3 min	20 min
World Health Organization (WHO)**	Pre-vacuum	134°C / 273°F	18 min	20 min

NB:

- * Steam sterilization cycles with longer times than those listed are also acceptable
- */** This cycle not recommended in the USA
- ** Recommended by the for reprocessing instruments where there is concern regarding TSE/CJD contamination. In the case of an operation with a risk of Non-Conventional Transmissible Agents (prions), this cycle must be completed by a period of total deactivation for the instruments concerned. Cf. Recommendations of the instruction DGS/RI3/2011/N°449 of 1st December 2011

Handling/Storage

- When handling of sterile boxes, special care must be paid so to avoid damaging the sterile barrier
- The health care establishment must establish a storage period for packaged instruments, depending on the sterile packaging and the recommendations of the manufacturer of this packaging
- The health care establishment must store sterile ancillaries in a suitable location which will allow damage to the sterile barrier to be prevented
- Instruments and ancillaries must not be exposed to direct sunlight, to ionising radiation or to extreme temperatures. They must be stored in a clean, ventilated room
- It is essential that a sterile ancillary can be distinguished from a non-sterile ancillary
- ADSM prohibits the transportation of sterile ancillaries outside the health care facility and disclaims any liability for such transportation.

References

- The instruction DGS/RI3/2011/N°449 of 1st December 2011 concerning the update of the recommendations aimed at reducing the risks of Non-Conventional Transmissible Agents (prions) during invasive operations
- **29CFR1910.1030**: Bloodborne pathogen
- **ISO 17664-1**: Sterilisation of medical devices -Information to be provided by the manufacturer for the processing of re-sterilisable medical devices
- **AAMI TIR12**: Designing, testing, and labelling reusable medical devices for re-processing in health care facilities: A guide for medical device manufacturers.
- **AAMI TIR 30**: A compendium of processes, materials, tests methods, and acceptance criteria for cleaning reusable devices
- **ANSI/AAMI ST79**: Comprehensive guide to steam sterilisation and sterility assurance in health care facilities.

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- **ANSI/AAMI ST35:** Safe handling and biological decontamination of reusable medical devices in health care facilities and in non-clinical settings
- **ISO 17665-1:** Sterilisation of health care products - Moist heat - Part 1: Requirements for the development, validation and routine control of a sterilisation process for medical devices
- **ISO 15883-1:** Washer-disinfectors -- Part 1: General requirements, terms and definitions and tests
- **ISO 15883-2:** Washer-disinfectors -- Part 2: Requirements and tests for washer-disinfectors employing thermal disinfection for surgical instruments, anaesthetic equipment, bowls, dishes, receivers, utensils, glassware, etc.

Contact

For any additional information, technical explanations, advice or assistance regarding this notice, contact:

7b rue Lavoisier
69680 Chassieu - FRANCE
Ph. : +33 (0)4 28 71 03 10
Fax : +33 (0)4 28 71 03 20
office@synchronomedical.com

**UK
CA**

CE₀₁₉₇

UK-RP : MDSS-UK RP Ltd.
6 Wilmslow Road, Rusholme
Manchester M14 5TP
United Kingdom

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Referenties

- The instruction **DGS/RI3/2011/N°449** of 1st December 2011 concerning the update of the recommendations aimed at reducing the risks of Non-Conventional Transmissible Agents (prions) during invasive operations
- **29CFR1910.1030**: Bloodborne pathogen
- **ISO 17664-1**: Sterilisation of medical devices -Information to be provided by the manufacturer for the processing of re-sterilisable medical devices
- **AAMI TIR12**: Designing, testing, and labelling reusable medical devices for re-processing in health care facilities: A guide for medical device manufacturers.
- **AAMI TIR 30**: A compendium of processes, materials, tests methods, and acceptance criteria for cleaning reusable devices
- **ANSI/AAMI ST79**: Comprehensive guide to steam sterilisation and sterility assurance in health care facilities.
- **ANSI/AAMI ST35**: Safe handling and biological decontamination of reusable medical devices in health care facilities and in non-clinical settings
- **ISO 17665-1**: Sterilisation of health care products - Moist heat - Part 1: Requirements for the development, validation and routine control of a sterilisation process for medical devices
- **ISO 15883-1**: Washer-disinfectors -- Part 1: General requirements, terms and definitions and tests
- **ISO 15883-2**: Washer-disinfectors -- Part 2: Requirements and tests for washer-disinfectors employing thermal disinfection for surgical instruments, anaesthetic equipment, bowls, dishes, receivers, utensils, glassware, etc.

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Viitteet

- The instruction DGS/RI3/2011/N°449 of 1st December 2011 concerning the update of the recommendations aimed at reducing the risks of Non-Conventional Transmissible Agents (prions) during invasive operations
- **29CFR1910.1030**: Bloodborne pathogen
- **ISO 17664-1**: Sterilisation of medical devices -Information to be provided by the manufacturer for the processing of re-sterilisable medical devices
- **AAMI TIR12**: Designing, testing, and labelling reusable medical devices for re-processing in health care facilities: A guide for medical device manufacturers.
- **AAMI TIR 30**: A compendium of processes, materials, tests methods, and acceptance criteria for cleaning reusable devices
- **ANSI/AAMI ST79**: Comprehensive guide to steam sterilisation and sterility assurance in health care facilities.
- **ANSI/AAMI ST35**: Safe handling and biological decontamination of reusable medical devices in health care facilities and in non-clinical settings
- **ISO 17665-1**: Sterilisation of health care products - Moist heat - Part 1: Requirements for the development, validation and routine control of a sterilisation process for medical devices
- **ISO 15883-1**: Washer-disinfectors -- Part 1: General requirements, terms and definitions and tests
- **ISO 15883-2**: Washer-disinfectors -- Part 2: Requirements and tests for washer-disinfectors employing thermal disinfection for surgical instruments, anaesthetic equipment, bowls, dishes, receivers, utensils, glassware, etc.

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