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INFORMATION MEDICALE IMPORTANTE TOEGRIP®CLASSIC et TOEGRIP®EVO

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INSTRUCTION FOR USE

IMPORTANT MEDICAL INFORMATION
TOEGRIP®CLASSIC and TOEGRIP®EVO

Reference: **SUP_7.015 version 4**
Date of last update: **05/2019**

Date of first EC marking:
- TOEGRIP®CLASSIC : 2010
- TOEGRIP®EVO : 2016

Product manufactured by:

ADSM
6 rue des Merisiers, 68920 Wettsolshem les Erlen – France
Tel: +33 (0)3 89 23 75 64 // Fax: +33 (0)9 70 19 01 32
adv@synchro-medical.com

CAUTION: USA Federal Law restricts this device to sale by or on the order of physician.

INTENDED USE – FDA Specific

The TOEGRIP®CLASSIC and TOEGRIP®EVO are intended for toes for general use in skeletally mature individuals undergoing surgery limited to interdigital fusion. The TOEGRIP®CLASSIC and TOEGRIP®EVO devices are intended to be permanently implanted without any other additional device as an intramedullary bone fastener device for **forefoot**.

PURPOSE

The TOEGRIP®CLASSIC and TOEGRIP®EVO interphalangeal implants are designed to relieve pain and disability of the forefoot by restoring and/or maintaining the alignment of two adjacent digital bone segments to optimize the achievement of a correct bone fusion for the concerned segments. These devices are only intended to be used for the forefoot of a mature skeleton.

DESCRIPTION:

The TOEGRIP®CLASSIC implant consists of a single part implant with three flexible intramedullary rods inserted into each part of the phalanx. The rods are attached by a T-shaped structure. The attachment concept is based on a press-fit contact due to its tapered shape and the macrostructures solidly anchored once impacted, thus avoiding any displacement of the implant. The TOEGRIP®CLASSIC implant is intended for single use only and is available in a range of 5 sizes with 3 possible degrees: 0°, 10° or 20°.

The TOEGRIP®EVO implant consists of a single part implant with three flexible intramedullary rods inserted into each part of the phalanx. The rods are connected by a cylindrical structure with a slight height offset between the proximal and distal portion. The attachment concept is based on a press-fit contact due to its tapered shape and the macrostructures solidly anchored once impacted, thus avoiding any displacement of the implant. The TOEGRIP®EVO implant is intended for single use only and is available in a range of 5 sizes with 3 possible degrees: 0°, 10° or 20°.

TOEGRIP®CLASSIC and TOEGRIP®EVO are made from PEEK according to Standard ASTM F2026.

INDICATIONS, CONTRAINDICATIONS AND POSSIBLE SIDE EFFECTS

Indications:
TOEGRIP®CLASSIC and TOEGRIP®EVO implants are indicated for the reconstruction of minor bones limited to the interdigital fusion of the toes on the **mature skeleton** in the following cases:

- **Rigidified deformations:**
 - Affected the proximal interphalangeal joint (PIJ).
 - Hammer toe deformities.
 - Claw toe deformities (PIP and DIP joints).
 - Hammer toe revision surgeries
 - Shortening osteotomies of the proximal phalanx

ATTENTION: to be used by or on the order of a licensed physician.

The licensed physician must take note of the documents accompanying the device. (IFU and surgical technique). No specific training is required for the understanding and implantation of the device. The medical doctor's qualifications and the reading of the accompanying documents are sufficient.

Contraindications:

- Non-exhaustive list of contraindications is as follows:
 - Any sign of generalised or local infection.
 - Pathological obesity.
 - Pregnancy.
 - Any other medical or surgical condition that may compromise the success of surgery with instruments, such as the presence of malignant tumours, or serious congenital anomalies, an increase in sedimentation rates that cannot be attributed to other diseases, an increase in the number of white blood cells or a downward trend in such blood cells.
 - Suspected or known allergy or intolerance to the implant's component materials.
 - Any situation requiring the use of different materials
 - Any case not listed in the indications.
 - Any patient who is not willing to follow the postoperative instructions.
 - Any patient in whom use of the implant may interfere with anatomical organs or some expected physiological function. The contraindications related to these devices are similar to those related to other interphalangeal devices. These interphalangeal instruments have not been designed for, intended or sold for any use other than those indicated.

Possible side effects:

- In addition to risks associated with any interphalangeal arthrodesis surgery, the list of possible undesirable side effects is as follows (non-exhaustive list):
- Migraion, fracture / rupture, expulsion, dislocation, or undesirable movement of one or more of the components
 - Persistent pain and disability
 - Bone loss related to short-circuiting of the constraints.
 - The body's reaction to a foreign object owing to the presence of the implant (e.g. the development of a tumour, auto-immune disease and / or poor healing)
 - Skin incompatibility with components where there is insufficient tissue coverage with pain and abnormal sensations due to the volume of the material.
 - Superficial or deep infection with inflammatory reaction.
 - Sepsis
 - Pseudarthrosis or delay of bone fusion.
 - Recurrence of the deformation, loss of correction
 - Excessive loss of blood during the operation and / or haematoma.
 - Lesions and / or vascular disorders (thrombosis)
 - Peripheral circulatory complication (thrombosis) and / or pulmonary embolism.
 - Inability to resume the activities of normal everyday life.
 - Compression of nearby tissues or organs.
- NOTE: Some additional surgical procedure may be required to remedy some of these possible undesirable events.

INSERTION OF THE DEVICE

The TOEGRIP®CLASSIC and TOEGRIP®EVO implants must be used with instruments specially designed for the device according to the surgical technique concerned. It must not be used with any other instruments, unless this has specifically been recommended in some other ADSM documentation, as using the system with other instruments poses the risk of incompatibility and cannot be guaranteed. ADSM shall not be held responsible in the case of using instruments that are not related to the implanted system. NEVER REUSE AN IMPLANT.
Reuse of an implant presents the risk of contamination and loss of mechanical integrity.

PACKAGING

The packaging of the device must be intact on receipt. Any device that has been damaged or any device whose packaging is damaged must not be used and must rather be returned to ADSM.

STORAGE AND MAINTENANCE CONDITIONS

The implants must be stored in their original packaging, out of direct sunlight, humidity and extreme temperatures.

MAGNETIC RESONANCE

The TOEGRIP®CLASSIC and TOEGRIP®EVO system has not been evaluated for safety and compatibility in the MR environment. The TOEGRIP®CLASSIC and TOEGRIP®EVO system has not been tested for heating or migration in the MR environment.

VERIFICATION

The devices must always be checked before use. Those that show signs of damage or scratches on the surface must not be used. The TOEGRIP®CLASSIC and TOEGRIP®EVO implants are delivered sterile and have been sterilised through gamma radiation. You need to carefully check that the round disc (sterilisation mark) on the packaging is red. Do not use an implant whose disc isn't red. You also need to check that the expiry date on the packaging has not passed. ADSM disengages itself from any liability in the case of use of implants past their expiry date. Any implant with damaged or torn packaging must not be resterilised but must rather be returned to ADSM.

TAKE NOTE: Do not resterilise an implant that has already been sterilised.

WARNINGS AND PRECAUTIONS

Every surgery carries the risk of failure, as many external factors can compromise the outcome, especially if the patient's state of health is already compromised. The implants must only be used in cases where bone fusion is deemed to be possible. No implant has been designed to withstand body weight in the long term without being accompanied by solid bone fusion. If a delay or lack of consolidation occurs, the implant may migrate and/or be damaged. **The functional product lifetime is of 6 months. In fact, the first bone fusion occurs around 90 days and keeping progress. Therefore, from 6 months the product is still implanted but it has no functional performance.** Given the current state of knowledge and device anteriority, **the risk caused by aging cannot be totally ruled out.**

The following of preoperative and operative procedures, knowledge of surgical techniques, the selection of implants and their correct placement as well as an appropriate selection and the patient's cooperation are factors leading to the operation's success. It has been shown that patients who smoke present a high risk of failure to obtain bone fusion. They need to be informed of this and any possible consequences. Obese patients, those who are malnourished and/or consume alcohol as well as patients with poor muscle or bone quality and/or those suffering from paralysis are also poor candidates for interphalangeal osteosynthesis surgery.

OTHER PREOPERATIVE, PEROPERATIVE AND POSTOPERATIVE WARNINGS

Implant selection
Selecting the suitable size of the implant for each patient is a crucial factor in the success of the operation, and it is important to refer to the surgical technique in this regard. Once implanted, the implant is subjected to repeated strain, and its strength is limited by adaptation of its geometrical structure to the size and shape of human bone. In order to minimise this strain and to avoid compromising the desired bone fusion, it is important to pay careful attention to the patient's selection criteria, the correct placement of the implant, and postoperative care.

In fact, conversely, these constraints can lead to excessive strain placed on the material, resulting in deformation, rupture or failure of the device, likely to cause damage or the need to remove the implant prematurely. Any use in body areas not recommended by ADSM cannot be guaranteed.

Preoperative precautions

- Only patients who meet the criteria described in the indications should be selected.
- Patients who meet the criteria described in the contraindications should not be selected.
- The implants must be handled and stored with the utmost care. They should not be scratched or damaged.
- Before carrying out the procedure, the surgeon must check availability of the range of implants and associated instruments. They must personally handle all the devices before use in order to familiarise themselves with insertion of the implant.
- Additional sterile components must be available in the event of unexpected need.

Preoperative precautions

- Instructions concerning the surgical technique should be followed carefully. Any breakage or improper use of the implant may injure the patient or surgical team.
- An implant holder must be used to insert the device.
- Take note: Do not use a device of the incorrect size (diameter and/or angulation), as this may damage the surrounding tissue.
- Where possible and necessary, an imaging system may be used to check on placement of the device.

Postoperative precautions

The instructions and postoperative warnings provided to the patient by the doctor, and the patient's ability to follow them are extremely important, so as not to compromise the desired bone fusion.

- Detailed instructions on the device's limits need to be provided to the patient. The patient must be informed that any deformity, loosening, and/or rupture of the device are complications that may occur during the postoperative rehabilitation phase following early use or excessive muscle activity. The risk of the onset of these complications could be heightened if the patient is active, frail, mentally ill or incapable of using some immobilisation device or external support. The patient must be advised in order to prevent falls, jumps or hard impacts on the area operated on or about any other movement that could compromise the desired bone fusion.
- To increase the likelihood of success of the surgery, the patient must not expose themselves to excessive strain which could cause weakening or rupture of the device. The patient must be informed of this risk and advised so they are aware of how to limit their physical activity, in particular twisting movements as well as any other participation in some sports activity.
- The patient should be informed about avoiding smoking or consuming alcohol during the bone consolidation process.
- Lack of bone fusion would lead to excessive, repeated strain placed on the implant. Through the mechanism of fatigue, this strain could end up causing some deformation, loosening or rupturing of the device. It is important to immobilise the fusion area and to check consolidation through x-ray. Should lack of bone fusion persist, or in the case of rupture, the device will need to be adjusted and/or removed immediately, before any serious injury occurs.
- An implant collected following ablation must be treated in a manner to prevent its reuse in some other surgical procedure, in accordance with the regulations described in the removal and disposal section for medical devices contained in this manual.

ADDITIONAL INFORMATION

For any additional information about the device or to enquire about specific surgical techniques, kindly contact ADSM's client services or the distributor.

PRODUCT COMPLAINTS

Any customer or user of this system who has a complaint or any reason for dissatisfaction concerning the quality of a product, its identity, durability, reliability, safety, effectiveness and/or its performance must notify the distributor or ADSM. In addition, should any of the implanted components not function well, or if poor functioning is suspected, the distributor or ADSM must be notified about this. If one of ADSM's products is ever suspected of having caused or contributed towards a patient's death or serious injury, ADSM or the distributor should be informed of this immediately by telephone, fax or written correspondence. For any complaint, please indicate the name, reference number and batch number of the component(s) concerned, your name and address, the nature of your complaint, and specify whether a written report is requested from the distributor or ADSM.

INFORMATION TO BE PROVIDED TO THE PATIENT

The TOEGRIP®CLASSIC and TOEGRIP®EVO interphalangeal implant is designed to restore and/or maintain two adjacent digital bone segments in anatomically correct alignment so as to optimise the obtaining of correct bone fusion for the segments concerned. This device is only intended for use in the forefoot. Specific instructions for the implant are described in the paragraph: «INDICATIONS, CONTRAINDICATIONS AND POSSIBLE UNDESIRABLE EFFECTS».

This device does not restore the functionality that can be expected of healthy bone and the patient should not expect unrealistic functional results. In addition, the anatomy of the human body limits the size of any artificial consolidation device used in surgery. This geometrical limitation increases the possibility of the onset of mechanical complications such as dismantling, deformation or rupturing of the device. Any complication may lead to an additional surgical procedure to remove the device or potentially insert another device. Consequently, it is very important that you carefully follow your doctor's postoperative instructions. It is recommended that the activities you do be limited to those on your surgeon's advice. Immobilisation and other devices that partially or fully support your weight must be used on your medical practitioner's advice. By following these instructions, you will increase the likelihood of achieving the desired result, and you will be reducing the risk of injury and/or additional surgical procedures.

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