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### INFORMATION MEDICALE IMPORTANTE

**TIGERTRACK<sup>®</sup> system**

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**CAUTION: USA Federal law restricts this device to sale by or on the order of physician.**

### INTENDED USE – FDA specific

The TIGERTRACK<sup>®</sup> screws system is intended for small bone extremely and largebone bone reconstruction fixation/arthrodesis procedures for general use in skeletally mature individuals. The TIGERTRACK<sup>®</sup> screws system is intended to be permanently implanted without any other additional device and are delivered sterile.

### PURPOSE

TIGERTRACK<sup>®</sup> screws are designed to relieve pain and disability of the forefoot by attaching and stabilising bone segments in elective osteotomies on metatarsal bones and phalanges of the foot, in order to optimise the correct bone fusion of the segments concerned. These devices are only intended for use in the forefoot of a mature skeleton.

### DESCRIPTION

The TIGERTRACK® screw system is intended for single use and is composed of TIGERTRACK® self-compressive headless screws (TIGER22, TIGER29 and TIGER45), TIGERTRACK® snap-off screws (TIGER18 and TIGER25) and Weil's snap-off screws TIGERW.

TIGERTRACK<sup>®</sup> self-compressive headless screws (TIGER22, TIGER29 and TIGER45) comprise a single part implant. These self-compressive headless screws are available in several different diameters and lengths.

Snap-off TIGERTRACK<sup>®</sup> screws (TIGER18 and TIGER25) comprise a single part implant fitted with a threaded portion and a snap-off connector portion with a wider diameter. These are self-compressive, headless screws. Weil's snap-off TIGERTRACK<sup>®</sup> screws (TIGERW) comprise a single part implant fitted with a threaded portion and a snap-off connector portion with a wider diameter. These are self-compressive screws with a head.

TIGERTRACK<sup>®</sup> screws are manufactured in Ti-6Al-4V ELI titanium according to Standard ISO 5832-3 / ASTM F 136. With the exception of the TIGER18 which is manufactured in a chrome-cobalt alloy according to ISO 5832-7.

### INDICATIONS, CONTRAINDICATIONS AND POSSIBLE SIDE EFFECTS

**Indications :**

The TIGERTRACK<sup>®</sup> screw system is indicated as an attachment system for small bone fractures or minor bone reconstruction in skeletally mature in the following cases:

- monocortical or bicortical osteotomy of the foot or hand
- metatarsal or distal or proximal metatarsal osteotomy
- fusion of the first metatarsophalangeal and interphalangeal joint
- attachment in the case of treatment of Hallux Valgus (e.g. scarf procedure, chevron osteotomy, etc.)
- Akin osteotomy
- Weil osteotomy

**CAUTION : to be used by or on the order of a licensed physician.**

The licensed physician must take note of the documents accompanying the device (IFU and surgical technique). No specific training is required for the understanding and implantation of the device. The medical doctor's qualifications and the reading of the accompanying documents are sufficient.

#### Contraindications:

A non-exhaustive list of contraindications is as follows:

1. Any sign of generalised or local infection.
  2. Pathological obesity.
  3. Pregnancy.
  4. Any other medical or surgical condition that may compromise the success of surgery with instruments, such as the presence of malignant tumours, or serious congenital anomalies, an increase in sedimentation rates that cannot be attributed to other diseases, an increase in the number of white blood cells or a downward trend in such blood cells.
  5. Suspected or known allergy or intolerance to the implant's component materials.
  6. Any situation requiring the use of different materials (including, but not listed in the indications).
  6. Any patient who is not willing to follow the postoperative instructions.
  9. Any patient in whom use of the implant may interfere with anatomical organs or some expected physiological function.
  10. Any dental treatment that could affect the passivation layer of metal implants (such as fluoride treatment)
- The contraindications pertaining to these devices are similar to those pertaining to other osteosynthesis instruments. These instruments have not been designed for, intended or sold for any use other than those indicated.

#### Possible side effects:

In addition to risks associated with any surgery to the extremities, the list of possible undesirable side effects is as follows (non-exhaustive list):

1. Migration, fracture/rupture, expulsion, dislocation, or undesirable movement of one or more of the components.
  2. Persistent pain and disability.
  3. Bone loss related to short-circuiting of the constraints.
  4. The body's reaction to a foreign object owing to the presence of the implant (e.g. the development of a tumour, auto-immune disease and/or poor healing).
  5. Skin incompatibility with components where there is insufficient tissue coverage with pain and abnormal sensations due to the volume of the material.
  6. Superficial or deep infection with inflammatory reaction.
  7. Sepsis.
  8. Pseudarthrosis or delay of bone fusion.
  9. Recurrence of the deformation, loss of correction.
  10. Excessive loss of blood during the operation and/or haematoma.
  11. Lesions and/or vascular disorders (thrombosis)
  12. Peripheral circulatory complication (thrombosis) and/or pulmonary embolism.
  13. Inability to resume the activities of normal everyday life.
  14. Compression of nearby tissues or organs.
- NOTE: Some additional surgical procedures may be required to remedy some of these possible undesirable events.

### INSERTION OF THE DEVICE

The TIGERTRACK<sup>®</sup> screw system must be used with instruments specially designed for the device according to the surgical technique concerned. It must not be used with any other instruments, unless this has specifically been recommended in some other ADSM documentation, as using the system with other instruments poses the risk of incompatibility, and cannot be guaranteed. ADSM shall not be held responsible in the case of using instruments that are not associated with the implanted system. NEVER REUSE AN IMPLANT. Reuse of an implant presents the risk of contamination and loss of mechanical integrity.

### PACKAGING

The packaging of the device must be intact on receipt. Any device that has been damaged or any device whose packaging is damaged must not be used, and must rather be returned to ADSM.

### STORAGE AND HANDLING CONDITIONS

The implants must be stored in their original packaging, out of direct sunlight, humidity and extreme temperatures.

### MAGNETIC RESONANCE

Magnetic Resonance Imaging (MRI) is a medical imaging technique that provides non-invasive views of the body interior using a strong magnetic field. Consequently, any ferromagnetic metallic foreign body is to be proscribed which can cause burns or injuries by displacement of the body under the effect of attraction of the magnetic field. Titanium Ti6Al4V ELI used for TIGER22, TIGER25, TIGER29, TIGER45 and TIGERW screws as well as Phynox used for TIGER18 screws are non-ferromagnetic metallic materials which pose no known risk for possible migration due to exposure to the magnetic resonance environment. It is the same for the heating which has been observed only on conductive loops which is not the case of the TIGERTRACK® system. A delay of 6-8 weeks is recommended after implantation before performing this type of examination. Some artifacts can be observed

### VERIFICATION

The devices must always be checked before use. Those that present signs of damage or scratches on the surface must not be used. TIGERTRACK<sup>®</sup> implants are delivered sterile, and have been sterilised through gamma radiation. You need to carefully check that the round disc (sterilisation mark) on the packaging is red. Do not use an implant whose disc isn't red.

You also need to check that the expiry date on the packaging has not passed. ADSM disengages itself from any liability in the case of the use of implants past their expiry date. Any implant with damaged or torn packaging must not be resterilised and must rather be returned to ADSM. TAKE NOTE: Do not resterilise an implant that has already been sterilised.

### WARNINGS AND PRECAUTIONS

Every surgery carries the risk of failure, as many external factors can compromise the outcome, especially if the patient's state of health is already compromised. The implants must only be used in cases where bone fusion is deemed to be possible. No implant has been designed to withstand body weight in the long term without being accompanied by solid bone fusion. If a delay or lack of consolidation occurs, the implant may migrate and/or rupture. The functional product lifetime is of 6 months. In fact, the first bone fusion occurs around 90 days and keeping progress. Therefore, from 6 months the product is still implanted but it has no functional performance.

Given the current state of knowledge and device anteriority, the risk caused by aging cannot be totally ruled out. The following of preoperative and operative procedures, knowledge of surgical techniques, the selection of implants and their correct placement as well as an appropriate selection and the patient's cooperation are factors leading to the operation's success. It has been shown that patients who smoke present a high risk of failure to obtain bone fusion. They need to be informed of this and any possible consequences. Obese patients, those who are malnourished and/or consume alcohol as well as patients with poor muscle or bone quality and/or those suffering from paralysis are also poor candidates for osteosynthesis surgery.

### OTHER PREOPERATIVE, PEROOPERATIVE AND POSTOPERATIVE WARNINGS

#### Implant selection

Selecting the suitable size of the implant for each patient is a crucial factor in the success of the operation, and it is important to refer to the surgical technique in this regard. Once implanted, the implant is subjected to repeated strain, and its strength is limited by adaptation of its geometrical structure to the size and shape of human bone. In order to minimise this strain and to avoid compromising the desired bone fusion, it is important to pay careful attention to the patient's selection criteria, the correct placement of the implant, and postoperative care.

In fact, conversely, these constraints can lead to excessive strain placed on the material, resulting in deformation, rupture or loosening of the device likely to cause damage or the need to remove the implant prematurely. Any use in body areas not recommended by ADSM cannot be guaranteed.

#### Preoperative precautions