5. Additional sterile components must be available in the event of unexpected need.

Per-operative precautions

1. Instructions concerning the surgical technique should be followed closely. Any breakage or improper use of the medical device may injure the patient or surgical team.
2. Take note: Do not use a device of the incorrect size (diameter and/or length), as this may damage the surrounding tissue.
3. Where possible and necessary, an imaging system may be used to check on placement of the device.

Post-operative precautions

The instructions and postoperative warnings provided to the patient by the doctor, and the patient's ability to follow them are extremely important, so as not to compromise the desired bone fusion.
1. Detailed instructions on the device's limits need to be provided to the patient. The patient must be informed that any dislodging, loosening and/or rotation of the device can complications that may occur during the postoperative rehabilitation phase following early use or excessive muscle activity. The risk of the onset of these complications could be heightened if the patient is active, frail, mentally ill or incapable of using some immobilisation device or external support. The patient must be advised in order to prevent falls, jumps or hard impacts on the area operated on or about any other action that could compromise the device's integrity.
2. To increase the likelihood of success of the surgery, the patient must not expose themselves to excessive strain which could cause weakening or rupture of the device. The patient must be informed of this risk and be advised so they are aware of how to limit their physical activity, in particular twisting movements as well as any other participation in some sports activities.
3. The patient should be informed about avoiding smoking or consuming alcohol during the bone consolidation process.
4. Lack of bone fusion would lead to excessive, repeated strain placed on the medical device. Through the mechanism of fatigue, this strain could end up causing some deformation, loosening or rupturing of the device. It is important to immobilise the fusion area and to check consolidation via x-ray. Should lack of bone fusion persist, or in the case of rupture, the device will need to be adjusted and/or removed immediately, before any serious injury occurs.
5. Any medical device collected following ablation must be treated in a manner to prevent its removal in some other surgical procedure, in accordance with the regulations described in the removal and disposal section for medical devices contained in this manual.

PACKAGING & TRANSPORT

Sterile delivered medical device

Medical devices manufactured by ADSM are sold sterile, using a gamma radiation sterilisation technique. The sterilisation method is indicated on the package. Items that have been sterilised through radiation have been exposed to at least 25 kGy of gamma radiation. If a device has been sterilised, it will be provided in a non-sterile state and must be sterilised before use (see Processing section).

Check the integrity of the packaging and labelling before opening. The packaging of the device must be intact on receipt. Any device that has been damaged or any device whose packaging is damaged must not be used and must be returned to ADSM.

Sterilisation is guaranteed on condition that the packaging is not damaged (torn wrap, a missing label, dubious packaging, etc.) and before the expiry date. Do not use a product whose packaging has been opened or is damaged outside of the operating room. The internal packaging must be handled under sterile conditions (persons/instruments).

Non-sterile delivered medical device

When a loan or deposit system is used, it must be carefully checked before use that no medical device is damaged or has signs of wear. If a medical device is damaged or has signs of wear, it must not be used and must be returned to ADSM or the distributor. Special attention must be paid to medical device with sharp and/or cutting components. It is imperative that they be placed in a sturdy container to prevent any risk of cutting or puncturing the container.

HANDLING

Sterilised medical devices must be handled in accordance with the appropriate aseptic procedures. It is necessary to protect medical devices from any accidental contact they may have with the product-specific medical devices to limit the risk of damaging the medical device surface. The care unit must carry out visual checking of the medical device before use to check that they are not damaged. Any medical device showing signs of damage must not be used. Additional information about these medical devices can be obtained from ADSM or the distributor.

PROCESSING (non-sterile delivered device) or REPROCESSING (reusable device for multiple patients)

It is incumbent to health care facilities to ensure that the re-treatment is done using appropriate equipment and materials and that the personnel responsible for the re-treatment have been correctly trained in order to obtain the desired result. The equipment and the processes must be validated and regularly monitored.

WARNINGS

Surgical instrumentation must be manipulated and/or used by trained, qualified persons having read the present instruction notice.

Surgical instrumentation must be used in premises compatible with standard aseptic conditions and with current practices in health care establishments.

The instruments and instruments are delivered clean but NOT STERILE.

The transfer of an ancillary between two health care establishments is forbidden without prior return to ADSTM.

Certain ancillaries are constituted by temporary invasive devices such as drilbits, reamers, miltips, taps and probes. These instruments must be used with precaution.

The equipment must only be used for the function for which they are dedicated in operating techniques. Any deterioration of the equipment can lead to risks of malfunction.

GENERALITIES

Before being used, the instruments, whether new or not, must be treated carefully according to the instructions.

The instruments boxes do not constitute a sterile barrier; a sterile packaging must be used to keep them in a sterile condition.

An instrument is delivered non-sterile in transport containers. A transfer sheet (which specifies the non-sterile state of the ancillary) is supplied with each ancillary for checking at reception and then before sterilisation.

Instruments or ancillaries generally have a long service life; however, poor handling or unsuitable protection can rapidly reduce it. Instruments which are no longer efficient because of long use, poor handling or inappropriate maintenance, must be returned to ADSTM. Any instrument intended to be dedicated to ADSTM and the distributor and the competent authority of the Member State in which the user and/or patient is established must be notified about this. If an instrument is damaged or has been damaged or caused or contributed towards a patient's death or serious injury, ADSTM or the distributor and the competent authority of the Member State in which the user and/or patient is established must be informed about this immediately by telephone, fax or written correspondence. For any complaint, please indicate the name, reference number and batch number of the concerned item, contact details, the nature of your complaint, and specify whether a written report is requested from the distributor or ADSTM. Any return of instruments to ADSTM or a distributor must be accompanied by a certificate guaranteeing that they have been cleaned, disinfected and sterilised beforehand.

INFORMATION TO BE PROVIDED TO THE PATIENT (only for long-term use devices)

STERILEBITZ & MISTERI devices for long-term use are designed for fixation and stabilization of bone segments for bone fusion. Specific indications are described in the paragraph: "INDICATIONS, CONTRAINDICATIONS AND POSSIBLE UNDESIRABLE EFFECTS".

This device does not restore the functionality that can be expected of healthy bone and the patient should not expect unrealistic functional results. In addition, the anatomy of the human body limits the size of any artificial consolidation device used in surgery. This geometrical limitation increases the possibility of the onset of mechanical complications such as dismantling, deformation or rupturing of the device. Any complication may lead to an additional surgical procedure to remove the device or potentially insert another device. Complications may also lead to secondary complications such as adhesions, infections, etc. It is recommended that the advices you do listen to those on your surgeon's advice. Immobilisation and other devices that partially or fully support your weight must be used on your medical practitioner's advice. By following these instructions, you will increase your chances of achieving the desired result, and you will be reducing the risk of injury and/or additional surgical procedures.

Only agents with a neutral pH must be used for cleaning aluminium devices. Contact with strong alkaline detergents or with products containing caustic soda, iodine or chlorine must be avoided as aluminium and stainless steel can be chemically attacked and the device damaged. The user must always refer to and comply with the instructions provided by the manufacturer.

The use of cobalt-chromium alloys is delivered non-sterile in transport containers. A transfer sheet (which specifies the non-sterile state of the ancillary) is supplied with each ancillary for checking at reception and then before sterilisation.

Any instrument intended to be dedicated to ADSTM and the distributor and the competent authority of the Member State in which the user and/or patient is established must be notified about this. If an instrument is damaged or has been damaged or caused or contributed towards a patient's death or serious injury, ADSTM or the distributor and the competent authority of the Member State in which the user and/or patient is established must be informed about this immediately by telephone, fax or written correspondence. For any complaint, please indicate the name, reference number and batch number of the concerned item, contact details, the nature of your complaint, and specify whether a written report is requested from the distributor or ADSTM. Any return of instruments to ADSTM or a distributor must be accompanied by a certificate guaranteeing that they have been cleaned, disinfected and sterilised beforehand.

CLEANING

Cautions

The following precautions must be observed by hospital personnel who work with contaminated or potentially contaminated medical devices. Pointed or cutting devices must be handled with great care.

An instrument is delivered non-sterile in transport containers. A transfer sheet (which specifies the non-sterile state of the ancillary) is supplied with each ancillary for checking at reception and then before sterilisation.

Instruments or ancillaries generally have a long service life; however, poor handling or unsuitable protection can rapidly reduce it. Instruments which are no longer efficient because of long use, poor handling or inappropriate maintenance, must be returned to ADSTM. Any instrument intended to be dedicated to ADSTM and the distributor and the competent authority of the Member State in which the user and/or patient is established must be notified about this. If an instrument is damaged or has been damaged or caused or contributed towards a patient's death or serious injury, ADSTM or the distributor and the competent authority of the Member State in which the user and/or patient is established must be informed about this immediately by telephone, fax or written correspondence. For any complaint, please indicate the name, reference number and batch number of the concerned item, contact details, the nature of your complaint, and specify whether a written report is requested from the distributor or ADSTM. Any return of instruments to ADSTM or a distributor must be accompanied by a certificate guaranteeing that they have been cleaned, disinfected and sterilised beforehand.

INSTRUCTIONS

When the device will be used

Remove excess liquid and organic tissues on the instruments with a lint-free disposable pad. Place the devices in a tray of distilled water and cover with a damp cloth.

Remove the plastic films and potential protective foam.

Visually check that all the instruments cited on the attached check-list are present.

The instruments must be cleaned within 30 minutes following their use to minimise the risk of cross-infection.

Conformity and transport

The instruments must be sent to the procurement service in closed or covered containers in order to avoid any new risk of contamination.

When this is applicable, the instruments with several components must be disassembled for efficient cleaning. Make sure not to lose small screws and components.

When this is applicable, articulated instruments must be open.

Antiseptics agents must be prepared observing the dilution and the temperature recommended by the manufacturer. Softened tap water can be used for preparing the cleaning agents.

The use of recommended temperatures is important for optimal performances of the cleaning agents.

Comment: new cleaning solutions must be prepared.

Procedure for automatic cleaning/disinfection

Sterika (cross-cleaning)

Completely immerse the instruments in an alkaline solution having fungicide, bactericides and virucide powers and leave them to soak for at least 15 minutes at room temperature. We recommend the use of the product NEODISHER®SEPTO PRECLEAN ZP ([C] = 5 mL / L) or the brand DR WEIGERT®. Use a nylon flexible brush to softly brush the device until all visible dirt has been removed. Pay particular attention to asperites, cavities, contact surfaces, connectors and other difficult to clean areas. Action the articulated parts in order to get to the difficult access areas.

Remove the device from the alkaline solution and rinse abundantly with lukewarm water of network for at least 5 minutes. Carefully and effectively rinse cavities, holes and any other difficult access areas.

Cleaning:

Place the instruments in the tray of a suitable washing machine/autoclave and treat using the standard cycle of the instrument washing machine/autoclave as per the instructions from the manufacturer of the washing machine/autoclave.

Washing with a detergent and disinfectant solution suitable for the washing machine (observe the concentrations and the temperature of use of the products).

We recommend the use of the product NEODISHER®SEPTO CLEAN ([C]=5 mL/L) DR WEIGERT® brand.

Disinfection:

Thermal disinfection assuring a value of A_0 greater than 3000, i.e. A minimum of 3 min at 93°C (or any other equivalent time/temperature pair)

Rinse with demineralized water. We recommend the use of MEDIKLAR® SPECIAL drying active ([C] = 0.6 mg / L) of the DR WEIGERT® brand.

Comment: the ancillary equipment is designed to operate at room temperature. It is recommended not to use the equipment directly after steam sterilisation and to wait for the temperature to return to ambient. Before the intervention, it is the responsibility of the personnel of the operating theatre to check the proper coupling of the ADSM instruments with the endotracheal tube and the oxygen cylinder in the operating theatre.

Carefully inspect each device to ensure the elimination of any visible contamination. If there is any contamination, repeat the cleaning/disinfection process.

Check the action of movable parts (e.g. hinges, connectors, sliding parts, springs, etc.) over their complete amplitude.

Check that there is no deformation of the instruments, in particular instruments which are connected to motors.

When the instruments are assembled, check that the devices fit together correctly with any components in contact.

Check that there is no excessive play between the assembled components.

PACKAGING

Medical packaging or steam sterilisation sachets available commercially, can be used for wrapping individual instruments. The wrapping must be prepared using a protective double wrapping technique or an equivalent method.

For sterilisation, trays and boxes with covers can also be placed in an approved sterilisation container with sealing lid. Follow the sterilisation container manufacturer's instructions for insertion and the replacement of sterilisation filters in sterilisation containers.

Instrument transport and storage boxes and trays must be used in the following conditions:

All devices must be properly arranged in order to ensure contact of the steam with all instrument surfaces. The instruments must not be stacked or placed in close contact use the arranging supports or silicone mat provided for this.

The user must check that the contents of the instrument box do not spill out, once the devices are arranged in the box.

Only devices manufactured and/or distributed by ADSM must be placed in ADSM trays.

STERILISATION

The decontamination step is compulsory before sterilising. It is important that sufficient cleaning is done before sterilisation without which, the sterilisation can be compromised.

After cleaning, the instruments must be rapidly and meticulously disassembled and placed in their respective trays/locations in order to avoid any possible damage or contamination. The instrument trays and supports must have been decontaminated and cleaned according to the same arrangement as those for the instruments.

The instruments must be sterilised in a dry environment for internal procedures for reassembly, inspection and packaging of instruments in order to ensure the sterilising penetration of the steam and appropriate drying. The arrangements to make for the protection of sharp or potentially injurious edges of instruments must also be recommended by the health care establishment.

Steam sterilisation heat is the recommended method for the ADSM set of orthopaedic instruments. Refer to the table below for the parameters of the sterilisation cycles validated by ADSM ensuring a sterility assurance level (SAL) of 10⁻⁶.

Sterilisation methods using ethylene oxide or plasma must not be used.

The recommendations of the manufacturer of the steriliser must always be followed. During the sterilisation of several instrument sets in a sterilisation cycle, ensure that the maximum load indicated by the manufacturer is not exceeded.

The set of instruments must be correctly prepared and packaged in the trays and/or boxes to allow the steam to penetrate and enter in direct contact with all surfaces.

Validated parameters for steam sterilization:

Type of cycle	Temperature (°C) / (°F)	Exposure Time	Drying Time
US	Pre-vacuum 132°C / 270°F	4 min	20 min
International*	Pre-vacuum 134°C / 273°F	3 min	20 min
World Health Organization (WHO)**	Pre-vacuum 134°C / 273°F	18 min	20 min

NB:

* Steam sterilization cycles with longer times than those listed are also acceptable

** This cycle not recommended in the USA

*** Recommended by the reprocessing instruments where there is concern regarding TSE/CD contamination. In the case of an operation with a risk of Non-Contaminable Transmissible Agents (prions), this cycle must be completed by a period of total deactivation for the instruments concerned. Cf. Recommendations of the instruction DGS/R13/2011/N°449 of 1st December 2011

HANDLING/STORAGE

When handling of sterile boxes, special care must be paid so to avoid damaging the sterile barrier.

The health care establishment must establish a storage period for packaged instruments, depending on the sterile packaging and the recommendations of the manufacturer of this packaging.

The health care establishment must store sterile ancillaries in a suitable location which will allow damage to the sterile barrier to be prevented.

Instruments and ancillaries must not be exposed to direct sunlight, to ionising radiation or to extreme temperatures. They must be stored in a clean, ventilated room

It is essential that a sterile ancillary can be distinguished from a non-sterile ancillary

ADSM prohibits the transportation of sterile ancillaries outside the health care facility and discards any liability for such transportation.

REFERRALS

The instruction DGS/R13/2011/N°449 of 1st December 2011 concerning the update of the operational procedures of the European Union concerning the sterilisation of medical devices.

ANSI/ANSI ST79: Comprehensive guide to steam sterilisation and sterility assurance in health care facilities

ANSI/ANSI ST35: Safe handling and biological decontamination of reusable medical devices in health care facilities and in non-clinical settings

ISO 17664-1: Sterilisation of medical devices - Information to be provided by the manufacturer for the processing of re-sterilisable medical devices

AAMI TIR12: Designing, testing, and labelling reusable medical devices for re-processing in health care facilities

The following are the main requirements of the instruction:

1. The instruction must be clear and concise.

2. The instruction must be legible and durable.

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AAMI TIR12: Designing, testing, and labeling reusable medical devices for re-processing in health care facilities: A guide for medical device manufacturers

AAMI TIR 30: A compendium of processes, materials, test methods and acceptance criteria for cleaning reusable devices

ANSI/AAMI ST79: Comprehensive guide to steam sterilization and sterility assurance in health care facilities

ANSI/AAMI ST35: Safe handling and biological decontamination of reusable medical devices in health care facilities and in nonclinical settings

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