



# TigerTrack<sup>®</sup>



Flange Screws



SynchroMedical

*The interphalangeal reference*

# TigerTrack®

## Flange Screw System

### DESCRIPTION

The TIGERTRACK® flange screw system is delivered in single sterile packed for single use and consists of the following products:

#### TIGER29 FLANGE *self-compressive headless screws* :



- Compressive and cannulated screw
- Self-drilling and self-tapping screw
- Flanged at 30°

#### TIGER35 FLANGE *non compressive headless screws* :



















- Non compressive, cannulated and fully threaded screw
- Self-drilling and self-tapping screw
- Flanged at 60°

### ORDERING INFORMATION

#### Implants

RANGE	DIAMETER	MIN LENGTH	MAX LENGTH
TIGER 29 F	∅ 2.9 mm	19	33
TIGER 35 F	∅ 3.5 mm	21	55

#### Instruments

	HV9010	AO HANDLE
	HVT7F	T7 CANNULATED SCREWDRIVER
	HV9023	T7 SCREWDRIVER
	HVT9F	T9 CANNULATED SCREWDRIVER
	HV9005	GAUGE LENGTH
	HV9035	GAUGE LENGTH
	HV9024	SHORT STEP DRILL
	HV9026	LONG STEP DRILL
	HV9025	CANNULATED SHORT STEP DRILL
	HV9027	CANNULATED LONG STEP DRILL
	HV93523C	CANNULATED SHORT STEP DRILL
	HV93523L	CANNULATED LONG STEP DRILL
	KW10100TR	KIRSCHNER WIRE WITH TROCAR TIP
	KW12150TR	KIRSCHNER WIRE WITH TROCAR TIP
	HV9015	INSTRUMENT TRAY
	HV9020	DEPTH GAUGE

TIGERTRACK® flange screws are manufactured in Ti-6Al-4V ELI titanium according to standards ISO 5832-3 / ASTM F 136..