



# TigerTrack<sup>®</sup>



Self-  
compressive  
Snap-off and  
Weil Screws



SynchroMedical

*The interphalangeal reference*

# TigerTrack®

## Screw System

### DESCRIPTION

The TIGERTRACK® screw system is delivered in single sterile packed for single use and consists of the following products:

#### TIGER22/29/35/45 self-compressive headless screws :



- Compressive and cannulated screw
- Self-drilling and self-tapping screw


#### TIGER18 & TIGER25 Snap-off threaded pins :



- Snap-off threaded pin Ø1,8 and Ø2,5 mm
- Power drill insertion



#### TIGERW Weil screws :



- Self drilling with cortical stop
- By tilting the driver to the side, the screw will break-off at the groove
- Finish insertion with hexalobe screwdriver

TIGERTRACK® screws are manufactured in Ti-6Al-4V ELI titanium according to standards ISO 5832-3 / ASTM F 136. With the exception of the TIGER18 which is manufactured in a chrome-cobalt alloy according to ISO 5832-7.

### ORDERING INFORMATION

#### Implants

RANGE	DIAMETER	MIN LENGTH	MAX LENGTH
TIGER 22	Ø 2.2 mm	11	45
TIGER 29	Ø 2.9 mm	11	45
TIGER 35	Ø 3.5 mm	23	61
TIGER 45	Ø 4.9 mm	25	61
TIGER 18	Ø 1.8 mm	13	45
TIGER 25	Ø 2.5 mm	13	45
TIGER W	Ø 2.2 mm	9	21

#### Instruments

	HV9010	AO HANDLE
	HV9003	T7 CANNULATED SCREWDRIVER
	HV9004	1,5MM HOLLOW HEXAGONAL SCREWDRIVER
	HV9023	T7 SCREWDRIVER
	HV9034	T9 CANNULATED SCREWDRIVER
	HV9005	GAUGE LENGTH
	HV9016	SCREW EXTRACTOR
	HV9035	GAUGE LENGTH
	HV9024	SHORT STEP DRILL
	HV9026	LONG STEP DRILL
	HV9025	CANNULATED SHORT STEP DRILL
	HV9027	CANNULATED LONG STEP DRILL
	HV9032	CANNULATED SHORT STEP DRILL
	HV9033	CANNULATED LONG STEP DRILL
	KW10100TR	KIRSCHNER WIRE WITH TROCER TIP
	KW08100TR	
	HV9015	INSTRUMENT TRAY
	HV9020	DEPTH GAUGE

All rights reserved © Pictures non-contractual